




Restorative Endodontic Excavators Set x 10

Description

Restorative Endodontic Excavators Set x 10:

Excavator 815-A, 1.2mm, Excavator 816-B, 1.4mm, Excavator 63/64, 1.5mm, Excavator 65/66, 1mm, Excavator 17W, 1.2/1.25mm, Excavator 18W, 1.5mm, Excavator 19W, 2mm, Excavator 125/126, 2.5mm, Excavator 127/128, 2mm, Excavator 129/130, 1.8 mm for removal of carious dentine. Manufactured from high grade surgical stainless steel.