



## Universal Implant 50 Piece Kit

### Description

Universal Implant 50 Piece Kit **(Only currently available on back order)**

This Kit Includes:

16pcs Conical Drills with Stoppers

1.8mm / 2.2mm – S1

1.8mm / 2.2mm – S3

1.8mm / 2.2mm – S5



1.8mm / 2.2mm – S7

2.2mm / 2.7mm – S1

2.2mm / 2.7mm – S3

2.2mm / 2.7mm – S5

2.2mm / 2.7mm – S7

2.3mm / 3.3mm – S1

2.3mm / 3.3mm – S3

2.3mm / 3.3mm – S5

2.3mm / 3.3mm – S7

2.8mm / 4.5mm – S1

2.8mm / 4.5mm – S3

2.8mm / 4.5mm – S5

2.8mm / 4.5mm – S7

\*S = Stopper Length

Bone Expander Screws

1.6mm / 2.6mm



1.9mm / 3.0mm

2.9mm / 3.4mm

2.3mm / 3.8mm

2.9mm / 4.3mm

7 Pcs Titanium Coated Drills

Drill Bit 2.0mm L 18.50

Drill Bit 2.8mm L 18.50

Drill Bit 3.65mm L 18.50

Drill Bit 4.2mm L 18.50

Drill Bit 4.5mm L 18.50

Drill Bit 4.8mm L 18.50

Drill Bit 5.2mm L 18.50

8 Pcs Trephine Drills

Int: 2.0mm, Ext: 3.0mm

Int: 3.0mm, Ext: 4.0mm

Int: 4.0mm, Ext: 5.0mm



Int: 5.0mm, Ext: 6.0mm

Int: 6.0mm, Ext: 7.0mm

Int: 7.0mm, Ext: 8.0mm

Int: 8.0mm, Ext: 9.0mm

Int: 9.0mm, Ext: 10.0mm

Compression Screw 2.8 / 4.0mm

Compression Screw 4.0mm

Driver 2.5mm

Compression Screw 2.0 / 2.8mm

Compression Screw 2.8 / 3.4mm

Tissue Punch 4.5mm

Tissue Punch 5.8mm

Pilot Drill 2.0mm

Drill Extender

Ratchet Driver 2.5mm

Hand Driver 2.5mm



Torque Wrench 6.35mm (10-40 N-cm)

Depth Gauge

Sterilization Cassette

Dental Conical Drill Kit Contents with Super Sharpness Backup  
Cutting Flutes

Each drill contains a stopper preventing over-drilling, reducing the potential for sinus or nerve complications during osteotomies. The system eliminates the need to carefully follow drill depth markings and for frequent pauses to examine osteotomy depth. Safely drill to the exact depth needed, with zero mistakes. The kit contains a separate drill for each diameter and length of implants. The drills are multiple use and in colours. This kit is cost-effective, user-friendly and uniquely compatible with any implant system.