



Prosthetic Screwdriver 8 Piece Set

Description

Prosthetic Screwdriver 8 Piece Set. Dental Implant Screw Set Abutment Tools Driver Holder Kit Drivers Instruments. Manufactured with French Stainless Steel Tips and Adonised Colour Handles, holder stainless steel.

This Kit Includes:

Green : Driver 1.0 ANK fit Antylos

Black : Driver 1.2 OSS fit 3I,SIC

Violet : Driver 1.27 DEN fit Bego, Zimmer, Astra

Yellow : Driver 1.4 ICX fit ICX

Red : Driver R155 NOB fit NOB

Blue : Driver R165 ITI fit ITI



Golden : Driver FRI 0.9 fit FRI

Orange : Driver ANT 1.16 fit ANT

Holder for 8 Screwdrivers