



Restorative Excavators Set x 10 – Hollow

Description

Restorative Endodontic Excavators Set x 10 Hollow Handles:
Excavator 1mm, Excavator 1.3mm Long Point, Excavator 1.3mm,
Excavator 1.5mm, Excavator 1.5mm Curved, Excavator 1.5mm
Curved Point, Excavator 1.8 mm, Excavator 2mm, Excavator 2.1mm,
Excavator 2.5mm for removal of carious dentine. Manufactured from
hollow high grade surgical stainless steel with French stainless steel
points.